



Apoptosis Detection Kits

PROMOTION

Apoptosis Detection Kits


25% Off

Limited-Time Promotion

Cell Apoptosis

Apoptosis: A Key Process in Programmed Cell Death, a form of programmed cell death (PCD), is a genetically controlled and orderly form of programmed cell death (PCD), is a genetically regulated process essential for cell growth, development, and renewal. Key events include phosphatidylserine (PS) externalization, mitochondrial membrane potential ($\Delta\Psi_m$) changes, caspase activation, and DNA fragmentation.


This orderly self-destruction occurs under physiological or pathological conditions, helping organisms adapt to their environment. Proper apoptosis is critical for maintaining cellular balance, while dysregulation can lead to diseases such as cancer and neurodegenerative disorders.



Elabscience® Apoptosis Detection Kits

Apoptosis detection is a critical tool for analyzing cell fate and disease mechanisms. Elabscience® offers comprehensive apoptosis detection kits covering multiple key apoptotic pathways, enabling precise identification of essential steps in the cell death process.

Annexin V Assay Kits




- 15 fluorochromes paired with 3 nuclear dyes
- Complete the experiment in just 15 minutes
- Distinct cell populations, effective even with high cell density

TUNEL Assay Kits




- Free of dimethyl arsenic, safe and non-toxic
- Wide variety, compatible with all detection equipment
- High sensitivity, comparable to leading international brands

Caspase Assay Kits



- Comprehensive panel, detects 8 types of Caspases
- Wide variety, available in both fluorescence and spectrophotometry methods
- Simple operation, delivering accurate experimental results

Mitochondrial Assay Kits




- Early apoptosis detection
- Simple experimental procedure
- Mitochondrial localization and function assessment



Apoptotic process events

1. The cell membrane loses its symmetry
2. Activation of apoptosis-related proteins such as Bcl-2
3. Caspase Activation
4. Mitochondrial membrane potential collapses, cytochrome C is released
5. Increase in sub-G1 cell population
6. Cell nuclear shrinkage
7. DNA fragmentation
8. Cell membrane budding to form apoptotic bodies



Schematic Diagram of Apoptosis Process Features



PROMOTION

Annexin V Assay Kits and Related Reagents

25% Off

PI

| Fluorochrome | Product Name |
|------------------------|---|
| FITC | Annexin V-FITC/PI Apoptosis Kit |
| Elab Fluor® 647 | Annexin V-Elab Fluor® 647/PI Apoptosis Kit |
| APC | Annexin V-APC/PI Apoptosis Kit |
| EGFP | Annexin V-EGFP/PI Apoptosis Kit |
| Cyanine5 | Annexin V-Cyanine5/PI Apoptosis Kit |
| APC/Cyanine7 | Annexin V-APC/Cyanine7/PI Apoptosis Kit |
| Elab Fluor® Violet 450 | Annexin V-Elab Fluor® Violet 450/PI Apoptosis Kit |
| Elab Fluor® Violet 500 | Annexin V-Elab Fluor® Violet 500/PI Apoptosis Kit |
| Elab Fluor® 488 | Annexin V-Elab Fluor® 488/PI Apoptosis Kit |
| Elab Fluor® Red 780 | Annexin V-Elab Fluor® Red 780/PI Apoptosis Kit |

7-AAD

| Fluorochrome | Product Name |
|------------------------|--|
| FITC | Annexin V-FITC/7-AAD Apoptosis Kit |
| Elab Fluor® 647 | Annexin V-Elab Fluor® 647/7-AAD Apoptosis Kit |
| PE | Annexin V-PE/7-AAD Apoptosis Kit |
| APC | Annexin V-APC/7-AAD Apoptosis Kit |
| EGFP | Annexin V-EGFP/7-AAD Apoptosis Kit |
| Cyanine5 | Annexin V-Cyanine5/7-AAD Apoptosis Kit |
| APC/Cyanine7 | Annexin V-PE/Cyanine7/7-AAD Apoptosis Kit |
| Elab Fluor® Violet 450 | Annexin V-Elab Fluor® Violet 450/7-AAD Apoptosis Kit |
| Elab Fluor® Violet 500 | Annexin V-Elab Fluor® Violet 500/7-AAD Apoptosis Kit |
| Elab Fluor® 488 | Annexin V-Elab Fluor® 488/7-AAD Apoptosis Kit |
| Elab Fluor® Red 780 | Annexin V-Elab Fluor® Red 780/7-AAD Apoptosis Kit |



DAPI

| Fluorochrome | Product Name |
|---------------------|--|
| FITC | Annexin V-FITC/DAPI Apoptosis Kit |
| Elab Fluor® 647 | Annexin V-Elab Fluor® 647/DAPI Apoptosis Kit |
| PE | Annexin V-PE/DAPI Apoptosis Kit |
| APC | Annexin V-APC/DAPI Apoptosis Kit |
| EGFP | Annexin V-EGFP/DAPI Apoptosis Kit |
| Cyanine5 | Annexin V-Cyanine5/DAPI Apoptosis Kit |
| PE/Cyanine5 | Annexin V-PE/Cyanine5/DAPI Apoptosis Kit |
| PE/Cyanine5.5 | Annexin V-PE/Cyanine5.5/DAPI Apoptosis Kit |
| PE/Cyanine7 | Annexin V-PE/Cyanine7/DAPI Apoptosis Kit |
| APC/Cyanine7 | Annexin V-APC/Cyanine7/DAPI Apoptosis Kit |
| PE/Elab Fluor® 594 | Annexin V-PE/Elab Fluor® 594/DAPI Apoptosis Kit |
| Elab Fluor® 488 | Annexin V-Elab Fluor® 488/DAPI Apoptosis Kit |
| Elab Fluor® Red 780 | Annexin V-Elab Fluor® Red 780/DAPI Apoptosis Kit |

Related Reagents

| Product Name |
|---|
| Recombinant Human ANXA5/Annexin V/Annexin A5 Protein(His tag) |
| ANXA5 Polyclonal Antibody |
| Annexin V Reagent |
| Annexin V Binding Buffer (10X) |
| Annexin V Lyophilized Powder |



PROMOTION

TUNEL Assay Kits

25% Off


Limited-Time Promotion

An important feature of the late stage of apoptosis is DNA fragmentation, and the exposed 3'-OH of the fragmented DNA can be catalyzed by terminal deoxynucleotidyl transferase (TdT) to bind to fluorochrome-labeled dUTP, which can be visualized by fluorescence microscopy or flow cytometry. This method is also known as One-step method (higher detection sensitivity, more intuitive results, and easier experimental results to be published in SCI articles).

If the exposed 3'-OH groups are labeled with biotin-labeled dUTP (Biotin-dUTP) under the catalysis of TdT, and then combined with HRP-labeled Streptavidin (Streptavidin-HRP), the result can be visualized by DAB color development under the catalysis of HRP, apoptotic cells can then be detected by an ordinary optical microscope.

Elabscience® developed the TUNEL series of apoptosis kits, including the **One-step TUNEL In Situ Apoptosis Kit, One-step TUNEL Flow Cytometry Apoptosis Kit, and TUNEL In Situ Apoptosis Kit (HRP-DAB Method)**. These kits can be used for apoptosis detection in tissue samples (paraffin-embedded, frozen sections) and cell samples (cell smears, cell culture films, suspension cells, adherent cells). Our TUNEL assay kits have high sensitivity, quick and easy operation, and can better assist in apoptosis research related to cell function and the R&D of related drugs.





One-step TUNEL In Situ Apoptosis Kits

Elabscience® One-step TUNEL In Situ Apoptosis Kits are suitable for in situ apoptosis detection of tissue samples (paraffin-embedded, frozen sections) and cell samples (cell smears, slide films), and the detection results can be directly observed by fluorescence microscopy. Compared with the DAB colorimetric method, the experimental results of fluorescent TUNEL are fluorescent images, which are more intuitive and easier to publish in SCI articles.



Features of Elabscience® One-step TUNEL In Situ Apoptosis Kits

Multiple Options

A variety of fluorescent dyes (FITC/ Elab Fluor® 488/ Elab Fluor® 594/ Elab Fluor® 647/ Elab Fluor® 555) kits are available.

Safe and Non-toxic

No dimethyl arsenate, non-toxic and harmless to people and the environment

Wide Range of Applications

Suitable for a variety of tissue samples (paraffin sections, frozen sections) and cell samples (cell slides)

Easy to Operate

One-step dyeing, saving time and effort.



Cost-effective

The key components are self-developed with more stable quality, and more cost-effective price

Elabscience® One-step TUNEL In Situ Apoptosis Kits

| Product Name | Dye | Excitation/emission wavelength | Size |
|---|-----------------------|---|------------------|
| One-step TUNEL In Situ Apoptosis Kit (Green, FITC) | FITC, DAPI | FITC: 495/525nm DAPI: 350/470nm | 20/50/100 Assays |
| One-step TUNEL In Situ Apoptosis Kit (Green, Elab Fluor® 488) | Elab Fluor® 488, DAPI | Elab Fluor® 488: 495/519nm DAPI: 350/470nm | 20/50/100 Assays |
| One-step TUNEL In Situ Apoptosis Kit (Red, Elab Fluor® 594) | Elab Fluor® 594, DAPI | Elab Fluor® 594: 590/617nm DAPI: 350/470nm | 20/50/100 Assays |
| One-step TUNEL In Situ Apoptosis Kit (Red, Elab Fluor® 647) | Elab Fluor® 647, DAPI | Elab Fluor® 647: 650/665nm DAPI: 350/470nm | 20/50/100 Assays |
| One-step TUNEL In Situ Apoptosis Kit (Red, Elab Fluor® 555) | Elab Fluor® 555, DAPI | Elab Fluor® 555: 555/565nm DAPI: 350/470nm | 20/100 Assays |

Comparative Results



Comparison and display of green fluorescence results between R company and Elabscience® Illumination value: 10, exposure time: 100 ms



One-step TUNEL Flow Cytometry Apoptosis Kits

Elabscience® One-step TUNEL Flow Cytometry Apoptosis Kit is a highly sensitive and quick and easy detection product for apoptosis. This kit is suitable for apoptosis detection of suspension cells and adherent cells, and the detection results can be analyzed by flow cytometry.

Features of Elabscience® One-step TUNEL Flow Cytometry Apoptosis Kits

Better detection scheme

Optimized fluorescent labeling solution, more suitable for flow cytometry detection schemes.

Easy to Operate

One-step dyeing, saving time and effort

Wide Range of Applications

Suitable for a variety of suspension cells and adherent cultured cells.

Multiple Options

FITC/Elab Fluor® 488/Elab Fluor® 594/Elab Fluor® Violet 450/Elab Fluor® 647/Elab Fluor® 555, a variety of fluorescence channels are available

Cost-effective


The kit contains fixed permeabilization fluid, no need to purchase additional reagents

Elabscience® One-step TUNEL Flow Cytometry Apoptosis Kit

| Product Name | Fluorochromes | Excitation/emission wavelength | Size |
|--|------------------------|--------------------------------|---------------|
| One-step TUNEL Flow Cytometry Apoptosis Kit (Green, FITC) | FITC | 495/525nm | 20/100 Assays |
| One-step TUNEL Flow Cytometry Apoptosis Kit (Green, Elab Fluor® 488) | Elab Fluor® 488 | 495/519nm | 20/100 Assays |
| One-step TUNEL Flow Cytometry Apoptosis Kit (Red, Elab Fluor® 594) | Elab Fluor® 594 | 590/617nm | 20/100 Assays |
| One-step TUNEL Flow Cytometry Apoptosis Kit (Blue, Elab Fluor® Violet 450) | Elab Fluor® Violet 450 | 410/450nm | 20/100 Assays |
| One-step TUNEL Flow Cytometry Apoptosis Kit (Red, Elab Fluor® 647) | Elab Fluor® 647 | 650/665nm | 20/100 Assays |
| One-step TUNEL Flow Cytometry Apoptosis Kit (Red, Elab Fluor® 555) | Elab Fluor® 555 | 555/565nm | 20/100 Assays |



Comparative Results



Biological control

Jurkat cells cultured without camptothecin, labeled with working solution, and TdT enzyme added

Negative control

Jurkat cells cultured with 5µM camptothecin for 4h, labeled with working solution, without TdT enzyme

Experimental group

Jurkat cells cultured with 5 µM camptothecin for 4h, labeled with working solution, and TdT enzyme added



TUNEL In Situ Apoptosis Kit (HRP-DAB Method)

Elabscience® TUNEL In Situ Apoptosis Kit (HRP-DAB Method) is a highly sensitive product that can detect apoptosis quickly and easily. This kit is suitable for tissue samples (paraffin sections, frozen sections). After biotin labeling and subsequent DAB color development, apoptosis can be observed with an ordinary optical microscope.

Features of Elabscience® TUNEL In Situ Apoptosis Kit (HRP-DAB Method)

High sensitivity

Very low background staining, strong positive staining and apoptosis can be detected at the single-cell level

Quick and easy

It only takes about 2~3 h to complete the staining

Wide range of applications

Suitable for a variety of tissue samples (paraffin sections, frozen sections)

Good specificity

It is generally easier to label apoptotic cells than necrotic cells in TUNEL assays

Safe and non-toxic


No dimethyl arsenate, non-toxic and harmless to people and the environment

Elabscience® TUNEL In Situ Apoptosis Kit (HRP-DAB Method)


| Cat.No. | Product Name | Size |
|-----------|--|------------------|
| E-CK-A331 | TUNEL In Situ Apoptosis Kit (HRP-DAB Method) | 20/50/100 Assays |



Results of TUNEL In Situ Apoptosis Kit (HRP-DAB Method)



Mouse Liver Tissue
Incubate with Dnase I for 10 min



Mouse Colon Tissue
Incubate with Dnase I for 10 min


PROMOTION

Caspase 25% Off

Limited-Time Promotion

Caspase, cysteinyl aspartate-specific proteinase, is a group of proteases present in the cytoplasm with similar structures. They exhibit high homology and structural similarity, capable of specifically cleaving peptide bonds at aspartate residues in target proteins. Caspases are closely associated with apoptosis (pyroptosis) in eukaryotic cells and participate in the regulation of cell growth, differentiation, and apoptosis.

Apoptosis



Elabscience®Caspase assay kits, for apoptosis/pyroptosis detection, high specificity, high sensitivity, simple and fast operation, excellent quality,can meet a wide range of your needs.




Detection

Colorimetric Method

Principle of Caspase Assay Kit (Colorimetric Method)

Caspase-specific peptide enzymes are linked to the yellow chromophore pNA (p-nitroanilide) as a substrate. Upon activation of Caspase, the enzyme cleaves the substrate, releasing the yellow chromophore. The absorbance of the released pNA can be measured using a spectrophotometer or enzyme reader at wavelengths around 405 nm or 400 nm, as pNA exhibits strong absorbance in this range.



Advantages

Multiple Choices

Provide a variety of Caspase assay kits to meet the different needs of customers

Cost-effective

The key components are developed in-house with reliable quality and more cost-effective price

Simple Operation

The composition of the reagent was optimized and the experimental operation was more time saving

Wide Applications

The enzyme activity in samples can be detected quantitatively, and these kits is suitable for a variety of cells

Caspase Assay Kits Colorimetric Method


| Product Name | Size | Application |
|--|---------------|---------------------|
| Caspase 1 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Cell Pyroptosis |
| Caspase 2 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Apoptosis Detection |
| Caspase 3/7 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Apoptosis Detection |
| Caspase 4 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Cell Pyroptosis |
| Caspase 6 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Apoptosis Detection |
| Caspase 8 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Apoptosis Detection |
| Caspase 9 Activity Assay Kit (Colorimetric Method) | 20/100 Assays | Apoptosis Detection |



Fluorescence Method

Principle of Caspase Assay Kit (Fluorescence Method)

Caspase-specific peptides are conjugated to high-affinity DNA fluorochrome as substrates. These substrates can penetrate the cell membrane and enter the cytoplasm. The substrates themselves are non-fluorescent and exhibit charge repulsion with DNA. Upon activation of Caspase, the enzyme cleaves the substrate, releasing the high-affinity DNA fluorochrome. This fluorochrome binds to DNA, resulting in strong fluorescence. This method allows for the detection of Caspase activity and observation of nuclear morphological changes during the process of apoptosis.



Advantages

Easy to operate

One-step dyeing, saving time and effort

Live-cell observation

Without lysis of the cells, live cell observations can be performed

Safe and non-toxic

Low cytotoxicity and did not affect the apoptotic process of cells

Wide Applications


Can be combined with multiple assay kits, choose diversification



Caspase Assay Kits (Fluorescence Method)

| Product Name | Cat.No. | Size |
|---|-----------|---------------|
| Caspase 1 Activity Detection Substrate for Flow Cytometry | E-CK-A481 | 20/100 Tests |
| Caspase 3/7 Activity Detection Substrate for Flow Cytometry | E-CK-A483 | 20/100 Tests |
| Caspase 4 Activity Detection Substrate for Flow Cytometry | E-CK-A484 | 20/100 Tests |
| Caspase 6 Activity Detection Substrate for Flow Cytometry | E-CK-A486 | 20/100 Tests |
| Caspase 8 Activity Detection Substrate for Flow Cytometry | E-CK-A488 | 20/100 Tests |
| Caspase 9 Activity Detection Substrate for Flow Cytometry | E-CK-A489 | 20/100 Tests |
| Caspase 3/7 and Annexin V Double Staining Apoptosis Kit | E-CK-A831 | 20/100 Assays |
| Caspase 3/7 and PI Double Staining Kit | E-CK-A832 | 20/100 Assays |
| Caspase 3/7 and DAPI Double Staining Kit | E-CK-A833 | 20/100 Assays |





www.microtech.eu

Sede Operativa
Traversa III Pisciarelli 18/a
80078 - Pozzuoli (NA)

Sede Legale
Viale Augusto 162,
80125 - Napoli (NA)

microtech@microtech.eu
microtech.srl@pec.it
+39 800 894 119

