




Puntali & Pipette

Puntali senza filtro

Compatibilità: Kartell®, Gilson®, Eppendorf®, Biohit®, Socorex®, HTL®, Nichiryo®, Finnpipette®, Thermo®

Realizzati in Polipropilene di elevata purezza ed esente da Cadmio. I puntali monouso vengono proposti in due differenti confezionamenti per soddisfare le diverse esigenze nella maniera maggiormente razionale:

- sfusi in sacchetto
- disposti in rack*

Scatole e puntali possono essere autoclavati a +121°C per 20 minuti.

*Le scatole che li contengono sono in Polipropilene e dotate di coperchio a pressione. Le scatole che contengono i puntali "macro" (5.000 µl e 10.000 µl) hanno il coperchio scorrevole orizzontalmente.

PUNTALI CAPACITÀ 0,1 - 10 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
956	Bulk	1000	52,50 €	20,00 €
9646	Rack (96)	960 (10x96)	170,00 €	61,50 €
PUNTALI CAPACITÀ 2-200 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
940	Bulk	1000	23,40 €	12,00 €
9643	Rack (96)	960 (10x96)	153,00 €	56,50 €
PUNTALI CAPACITÀ 50-1.000 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
941	Bulk	1000	32,10 €	12,00 €
9645	Rack (96)	960 (10x96)	153,00 €	59,50 €
PUNTALI CAPACITÀ 1.000-5.000 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
984	Bulk	250	57,80 €	57,80 €
9649	Rack (96)	200 (4x50)	169,20 €	169,20 €
PUNTALI CAPACITÀ 1.000-5.000 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
984	Bulk	250	57,80 €	57,80 €
9649	Rack (96)	200 (4x50)	169,20 €	169,20 €
PUNTALI CAPACITÀ 2.000-10.000 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
943	Bulk	100	52,00 €	52,00 €
9848	Rack (96)	200 (4x50)	189,20 €	189,20 €

Puntali con filtro

Compatibilità: Gilson®, Eppendorf®, Socorex®, HTL®, Nichiryo®, Finnpipte®, Thermo®

Puntali sterilizzati, in rack, realizzati in Polipropilene e dotati di filtro.
Sono totalmente esenti da DNase, RNase, DNA, ATP, metalli e Cadmio.
Apirogeni ed esenti da inibitori della PCR.

Sterilizzati con radiazioni Gamma, possono essere autoclavati a +121°C per 20 minuti.

Ideali per impieghi in biologia molecolare, colture cellulari, batteriologia, virologia, immunologia e radioimmunologia.

Item	Capacity (µl)	Size	List Price	Promo Price
9638	0,1-10	960 (10x96)	313,00 €	79,00 €
9634	02-30	960 (10x96)	313,00 €	79,00 €
9635	2-100	960 (10x96)	313,00 €	79,00 €
9629	2-200	960 (10x96)	313,00 €	79,00 €
9633	50-1000	960 (10x96)	313,00 €	79,00 €



Hipette single-channel adjustable volume pipettes

Introducing the new generation of fully autoclavable manual pipettes. The DLAB pipette.

- **Ergonomic Design:** The DLAB pipette is designed with excellent ergonomics that reduce fatigue and provide a smoother and easier pipetting experience.
- **Smooth and Effortless Operation:** With a one-stop point of 4N and two-stop points of 10N, the DLAB pipette offers a comfortable single-handed pipetting experience with smooth and effortless operation.
- **Multiple Volume Range:** The DLAB pipette offers 12 different single-channel adjustable volume , ranging from 0.1 µl to 10 ml, making it ideal for various laboratory applications.
- **Precise Volume Setting:** The DLAB pipette features a new piston and metal screw structure, ensuring precise suction and dispensing of liquids. Its counter increment is one thousandth of the volume range, making it easy to set and read.
- **One-Click Volume Lock:** The DLAB pipette has a one-click volume lock feature that helps prevent pipetting errors caused by accidental volume changes, allowing for more accurate transfer of small volumes of liquids.
- **Smooth Volume Adjustment:** The DLAB pipette offers a smooth and easy volume adjustment experience, ensuring high precision and accuracy.
- **High-Quality Materials:** The DLAB pipette is made of a new type of material that does not require disassembly for autoclaving or UV disinfection, and is resistant to harsh chemicals.

Item	Capacity (µl)	Packaging	Promo Price
7030511001	0.1 - 2.5	1	91,50 €
7030511004	0.5 - 10	1	91,50 €
7030511005	02 - 20	1	91,50 €
7030511008	10 - 100	1	91,50 €
7030511009	20 - 200	1	91,50 €
7030511040	30 - 300	1	91,50 €
7030511014	100 - 1000	1	91,50 €
7030301017	1000 - 5000	1	91,50 €



Gilson original not Filtered Pipette Tips

PIPETMAN® DIAMOND Tips are designed to fit your PIPETMAN pipettes perfectly. They deliver accurate and precise liquid volumes when used with Gilson pipettes.

EASYPACK™ is an economical way to order tips in bulk.

PUNTALI CAPACITÀ 0,1 - 10 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
F161631	Bulk EASYPACK	1000	49,00 €	40,00 €
F161630	Bulk ECOPACK	10X1000	393,00 €	315,00 €
F171100	Rack (96)	960 (10x96)	78,00 €	63,00 €
F171101	Rack (96) Sterile	960 (10x96)	87,00 €	70,00 €
F167101	Refil (96)	960 (10x96)	65,00 €	53,00 €
F167201	Refil (96) Sterile	960 (10x96)	80,00 €	65,00 €
PUNTALI CAPACITÀ 2-200 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
F161931	Bulk EASYPACK	1000	32,00 €	26,00 €
F161930	Bulk ECOPACK	10X1000	258,00 €	207,00 €
F171300	Rack (96)	960 (10x96)	65,00 €	53,00 €
F171301	Rack (96) Sterile	960 (10x96)	73,00 €	59,00 €
F167103	Refil (96)	960 (10x96)	58,00 €	47,00 €
F167203	Refil (96) Sterile	960 (10x96)	65,00 €	53,00 €
PUNTALI CAPACITÀ 50-1.000 microlitri				
Item	Packaging	Pack Size	List Price	Promo Price
F161671	Bulk ECOPACK	1000	32,00 €	26,00 €
F161670	Bulk EASYPACK	10X1000	258,00 €	207,00 €
F171500	Rack (96)	960 (10x96)	67,00 €	54,00 €
F171501	Rack (96) Sterile	960 (10x96)	78,00 €	63,00 €
F167104	Refil (96)	960 (10x96)	48,00 €	39,00 €
F167204	Refil (96) Sterile	960 (10x96)	60,00 €	49,00 €

Gilson original Sterilized Filter Tips

PIPETMAN® DIAMOND Sterilized Filter Tips are designed to fit your PIPETMAN perfectly and play a critical role in preventing aerosol-born contamination.

They deliver accurate and precise pipetting when used in conjunction with Gilson pipette. Ready to use - no need to autoclave, and color-coded for easy identification. Compliant with Good Laboratory Practice (GLP).

Item	Capacity (µl)	Pack Size	List Price	Promo Price
F171103	0,1-10	960 (10x96)	172,00 €	138,00 €
F171303	02-30	960 (10x96)	165,00 €	133,00 €
F171403	10-100	960 (10x96)	165,00 €	133,00 €
F171503	20-200	960 (10x96)	135,00 €	109,00 €
F171603	20-300	960 (10x96)	174,00 €	140,00 €
F171703	100-1000	960 (10x96)	194,00 €	156,00 €

Gilson PIPETMAN®, recognized as the pipetting standard for its reliability and proven durability, is a fully adjustable, air-displacement pipette that uses disposable tips. PIPETMAN combines legendary accuracy, precision, and durability with easy pipetting and handling. PIPETMAN pipettes come in a large range of volumes with reduced pipetting forces for hours of comfortable pipetting.

Item	Capacity (µl)	Model	List Price	Promo Price
F144054M	0,2-2	P2	323,00 €	259,00 €
F144055M	01-10	P10	323,00 €	259,00 €
F144056M	02-20	P20	323,00 €	259,00 €
F144057M	10-100	P100	323,00 €	259,00 €
F144058M	20-200	P200	323,00 €	259,00 €
F144058M	100-1000	P1000	323,00 €	259,00 €
F144066	1000-5000	P5000	323,00 €	259,00 €
F144067	1000-10000	P10mL	323,00 €	259,00 €



www.microtech.eu

Sede Operativa
Traversa III Pisciarelli 18/a
80078 - Pozzuoli (NA)

Sede Legale
Viale Augusto 162,
80125 - Napoli (NA)

microtech@microtech.eu
microtech.srl@pec.it
+39 800 894 119

