



# Puntali & Pipette



# Puntali senza filtro

**Compatibilità: Kartell®, Gilson®, Eppendorf®, Biohit® Socorex®, HTL®, Nichiryo®, Finnpiquette®, Thermo®**

Realizzati in Polipropilene di elevata purezza ed esente da Cadmio. I puntali monouso vengono proposti in due differenti confezionamenti per soddisfare le diverse esigenze nella maniera maggiormente razionale:

- sfusi in sacchetto
- disposti in rack\*

Scatole e puntali possono essere autoclavati a +121°C per 20 minuti.

\*Le scatole che li contengono sono in Polipropilene e dotate di coperchio a pressione. Le scatole che contengono i puntali "macro" (5.000 µl e 10.000 µl) hanno il coperchio scorrevole orizzontalmente.

| PUNTALI CAPACITÀ 0,1 - 10 microlitri     |           |             |            |             |
|------------------------------------------|-----------|-------------|------------|-------------|
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 956                                      | Bulk      | 1000        | 52,50 €    | 20,00 €     |
| 9646                                     | Rack (96) | 960 (10x96) | 170,00 €   | 61,50 €     |
| PUNTALI CAPACITÀ 2-200 microlitri        |           |             |            |             |
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 940                                      | Bulk      | 1000        | 23,40 €    | 12,00 €     |
| 9643                                     | Rack (96) | 960 (10x96) | 153,00 €   | 56,50 €     |
| PUNTALI CAPACITÀ 50-1.000 microlitri     |           |             |            |             |
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 941                                      | Bulk      | 1000        | 32,10 €    | 12,00 €     |
| 9645                                     | Rack (96) | 960 (10x96) | 153,00 €   | 59,50 €     |
| PUNTALI CAPACITÀ 1.000-5.000 microlitri  |           |             |            |             |
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 984                                      | Bulk      | 250         | 57,80 €    | 57,80 €     |
| 9649                                     | Rack (96) | 200 (4x50)  | 169,20 €   | 169,20 €    |
| PUNTALI CAPACITÀ 1.000-5.000 microlitri  |           |             |            |             |
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 984                                      | Bulk      | 250         | 57,80 €    | 57,80 €     |
| 9649                                     | Rack (96) | 200 (4x50)  | 169,20 €   | 169,20 €    |
| PUNTALI CAPACITÀ 2.000-10.000 microlitri |           |             |            |             |
| Item                                     | Packaging | Pack Size   | List Price | Promo Price |
| 943                                      | Bulk      | 100         | 52,00 €    | 52,00 €     |
| 9848                                     | Rack (96) | 200 (4x50)  | 189,20 €   | 189,20 €    |

# Puntali con filtro

**Compatibilità: Gilson®, Eppendorf®, Socorex®, HTL®, Nichiryo®, Finnpipette®, Thermo®**

Puntali sterilizzati, in rack, realizzati in Polipropilene e dotati di filtro.  
Sono totalmente esenti da DNase, RNase, DNA, ATP, metalli e Cadmio.  
Apirogeni ed esenti da inibitori della PCR.

Sterilizzati con radiazioni Gamma, possono essere autoclavati a +121°C per 20 minuti.

**Ideali per impieghi in biologia molecolare, colture cellulari, batteriologia, virologia, immunologia e radioimmunologia.**

| Item | Capacity (µl) | Size        | List Price | Promo Price |
|------|---------------|-------------|------------|-------------|
| 9638 | 0,1-10        | 960 (10x96) | 313,00 €   | 79,00 €     |
| 9634 | 02-30         | 960 (10x96) | 313,00 €   | 79,00 €     |
| 9635 | 2-100         | 960 (10x96) | 313,00 €   | 79,00 €     |
| 9629 | 2-200         | 960 (10x96) | 313,00 €   | 79,00 €     |
| 9633 | 50-1000       | 960 (10x96) | 313,00 €   | 79,00 €     |



# Hipette single-channel adjustable volume pipettes

Introducing the new generation of fully autoclavable manual pipettes. The DLAB pipette.

- **Ergonomic Design:** The DLAB pipette is designed with excellent ergonomics that reduce fatigue and provide a smoother and easier pipetting experience.
- **Smooth and Effortless Operation:** With a one-stop point of 4N and two-stop points of 10N, the DLAB pipette offers a comfortable single-handed pipetting experience with smooth and effortless operation.
- **Multiple Volume Range:** The DLAB pipette offers 12 different single-channel adjustable volume, ranging from 0.1  $\mu$ l to 10 ml, making it ideal for various laboratory applications.
- **Precise Volume Setting:** The DLAB pipette features a new piston and metal screw structure, ensuring precise suction and dispensing of liquids. Its counter increment is one thousandth of the volume range, making it easy to set and read.
- **One-Click Volume Lock:** The DLAB pipette has a one-click volume lock feature that helps prevent pipetting errors caused by accidental volume changes, allowing for more accurate transfer of small volumes of liquids.
- **Smooth Volume Adjustment:** The DLAB pipette offers a smooth and easy volume adjustment experience, ensuring high precision and accuracy.
- **High-Quality Materials:** The DLAB pipette is made of a new type of material that does not require disassembly for autoclaving or UV disinfection, and is resistant to harsh chemicals.

| Item       | Capacity ( $\mu$ l) | Packaging | Promo Price |
|------------|---------------------|-----------|-------------|
| 7030511001 | 0.1 - 2.5           | 1         | 91,50 €     |
| 7030511004 | 0.5 - 10            | 1         | 91,50 €     |
| 7030511005 | feb-20              | 1         | 91,50 €     |
| 7030511008 | 10 - 100            | 1         | 91,50 €     |
| 7030511009 | 20 - 200            | 1         | 91,50 €     |
| 7030511040 | 30 - 300            | 1         | 91,50 €     |
| 7030511014 | 100 - 1000          | 1         | 91,50 €     |
| 7030301017 | 1000 - 5000         | 1         | 91,50 €     |

# Gilson original not Filtered Pipette Tips

**PIPETMAN® DIAMOND Tips** are designed to fit your PIPETMAN pipettes perfectly. They deliver accurate and precise liquid volumes when used with Gilson pipettes.

EASYPACK™ is an economical way to order tips in bulk.

| PUNTALI CAPACITÀ 0,1 - 10 microlitri |                    |             |            |             |
|--------------------------------------|--------------------|-------------|------------|-------------|
| Item                                 | Packaging          | Pack Size   | List Price | Promo Price |
| F161631                              | Bulk EASYPACK      | 1000        | 49,00 €    | 40,00 €     |
| F161630                              | Bulk ECOPACK       | 10X1000     | 393,00 €   | 315,00 €    |
| F171100                              | Rack (96)          | 960 (10x96) | 78,00 €    | 63,00 €     |
| F171101                              | Rack (96) Sterile  | 960 (10x96) | 87,00 €    | 70,00 €     |
| F167101                              | Refil (96)         | 960 (10x96) | 65,00 €    | 53,00 €     |
| F167201                              | Refil (96) Sterile | 960 (10x96) | 80,00 €    | 65,00 €     |
| PUNTALI CAPACITÀ 2-200 microlitri    |                    |             |            |             |
| Item                                 | Packaging          | Pack Size   | List Price | Promo Price |
| F161931                              | Bulk EASYPACK      | 1000        | 32,00 €    | 26,00 €     |
| F161930                              | Bulk ECOPACK       | 10X1000     | 258,00 €   | 207,00 €    |
| F171300                              | Rack (96)          | 960 (10x96) | 65,00 €    | 53,00 €     |
| F171301                              | Rack (96) Sterile  | 960 (10x96) | 73,00 €    | 59,00 €     |
| F167103                              | Refil (96)         | 960 (10x96) | 58,00 €    | 47,00 €     |
| F167203                              | Refil (96) Sterile | 960 (10x96) | 65,00 €    | 53,00 €     |
| PUNTALI CAPACITÀ 50-1.000 microlitri |                    |             |            |             |
| Item                                 | Packaging          | Pack Size   | List Price | Promo Price |
| F161671                              | Bulk ECOPACK       | 1000        | 32,00 €    | 26,00 €     |
| F161670                              | Bulk EASYPACK      | 10X1000     | 258,00 €   | 207,00 €    |
| F171500                              | Rack (96)          | 960 (10x96) | 67,00 €    | 54,00 €     |
| F171501                              | Rack (96) Sterile  | 960 (10x96) | 78,00 €    | 63,00 €     |
| F167104                              | Refil (96)         | 960 (10x96) | 48,00 €    | 39,00 €     |
| F167204                              | Refil (96) Sterile | 960 (10x96) | 60,00 €    | 49,00 €     |



# Gilson original Sterilized Filter Tips

**PIPETMAN® DIAMOND Sterilized Filter Tips** are designed to fit your PIPETMAN perfectly and play a critical role in preventing aerosol-born contamination.

They deliver accurate and precise pipetting when used in conjunction with Gilson pipette. Ready to use - no need to autoclave, and color-coded for easy identification. Compliant with Good Laboratory Practice (GLP).

| Item    | Capacity (µl) | Pack Size   | List Price | Promo Price |
|---------|---------------|-------------|------------|-------------|
| F171103 | 0,1-10        | 960 (10x96) | 172,00 €   | 138,00 €    |
| F171303 | feb-30        | 960 (10x96) | 165,00 €   | 133,00 €    |
| F171403 | 10-100        | 960 (10x96) | 165,00 €   | 133,00 €    |
| F171503 | 20-200        | 960 (10x96) | 135,00 €   | 109,00 €    |
| F171603 | 20-300        | 960 (10x96) | 174,00 €   | 140,00 €    |
| F171703 | 100-1000      | 960 (10x96) | 194,00 €   | 156,00 €    |

Gilson PIPETMAN®, recognized as the pipetting standard for its reliability and proven durability, is a fully adjustable, air-displacement pipette that uses disposable tips. PIPETMAN combines legendary accuracy, precision, and durability with easy pipetting and handling. PIPETMAN pipettes come in a large range of volumes with reduced pipetting forces for hours of comfortable pipetting.

| Item     | Capacity (µl) | Model | List Price | Promo Price |
|----------|---------------|-------|------------|-------------|
| F144054M | 0,2-2         | P2    | 323,00 €   | 259,00 €    |
| F144055M | 01-10         | P10   | 323,00 €   | 259,00 €    |
| F144056M | 02-20         | P20   | 323,00 €   | 259,00 €    |
| F144057M | 10-100        | P100  | 323,00 €   | 259,00 €    |
| F144058M | 20-200        | P200  | 323,00 €   | 259,00 €    |
| F144058M | 100-1000      | P1000 | 323,00 €   | 259,00 €    |
| F144066  | 1000-5000     | P5000 | 323,00 €   | 259,00 €    |
| F144067  | 1000-10000    | P10mL | 323,00 €   | 259,00 €    |

**www.microtech.eu**

**Sede Operativa**

Traversa III Pisciarelli 18/A  
80078 - Pozzuoli (NA)

**Sede Legale**

Viale Augusto 162,  
80125 - Napoli (NA)

microtech@microtech.eu  
microtech.srl@pec.it  
+39 800 894 119

